

March 2016

ENTERAL CONNECTIONS: NEW STANDARD & DESIGNS

Medtronic
Further, Together

CONTENT

- ISO Standard 80369-3
- What is changing?
 - Feeding sets
 - Syringes
 - Feeding tubes
 - Packaging
- ANZ timelines
- Industry collaboration



KANGAROO CONNECTIONS

ISO STANDARD 80369-3

What is it?

- New ISO standards have been established for small bore connectors, not only on the proximal/nutritional end but also on the distal/patient access end.


What you can expect?

- Industry wide changes to the distal connection on all pump sets and a corresponding change to all feeding tube ports.

KANGAROO CONNECTIONS

WHAT IS CHANGING? FEEDING SETS

Present feeding sets



KANGAROO CONNECTIONS

WHAT IS CHANGING? FEEDING SETS

Transition feeding sets: 2015 to 2016

Silicone tubing to be loaded into the pump

Patient connector

The new ENFit stepped connector will be pre-attached to feeding set tubing during the transition period

Medication port

The new ENFit medication connector will be included loose in the packaging during the transition period

Pump



KANGAROO CONNECTIONS

WHAT IS CHANGING? FEEDING TUBES

Present feeding tubes



KANGAROO CONNECTIONS

WHAT IS CHANGING? FEEDING TUBES

Future feeding tubes: 2016 to 2017



KANGAROO CONNECTIONS

WHAT IS CHANGING? SYRINGES


Present enteral syringes



KANGAROO CONNECTIONS

WHAT IS CHANGING? SYRINGES

Future enteral syringes: 2017






Note: Medtronic will not be providing ENFit syringes in the short term. However, we will be providing syringe port adaptors, enabling the use of luer lock and luer slip syringes in medication ports. Please contact your existing syringe supplier for further information on ENFit syringes.

KANGAROO CONNECTIONS

WHAT IS CHANGING? PACKAGING

External packaging

Current		ENFit
Original Shipping Label 		ENFit Shipping Label 

Internal packaging





- Feeding sets that do not include medication ports will have yellow leaflets within the internal packaging
- Feeding sets that do include medication ports will remain unchanged, with white leaflets within the internal packaging.



KANGAROO CONNECTIONS

ANZ TIMELINES

ANZ Timelines

Timing		ENFit
2015 - 2016		Feeding Sets
2016		Transition Connectors
2016 - 2017		Feeding Tubes
2017		Syringes

THANK YOU
Q&A